

Snowden-Pencer[®] articulating instruments




Snowden-Pencer® articulating instruments


Snowden-Pencer articulating instruments use flexible technology that allows the instruments to act like a natural extension of a surgeon's hand.

The articulating graspers and dissectors are designed to transmit grasping force through the delicate working tips during laparoscopic and single port procedures. These instruments allow the surgeon to "reach around the corner" or dissect around difficult anatomy during demanding procedures.


Articulating forceps are available in both 40° and 90° configurations in five different grasper jaw patterns as well as a Maryland dissector jaw pattern.




Scissor-style handle




Articulating BABCOCK grasper, 5mm

	Catalog no.	Dimensions
	89-0504	40°, 34cm length
	89-0500	90°, 34cm length
	89-4504	40°, 45cm length
	89-4500	90°, 45cm length

Articulating DeBAKEY grasper, 5mm

	Catalog no.	Dimensions
	89-0506	40°, 34cm length
	89-0507	90°, 34cm length
	89-4506	40°, 45cm length
	89-4507	90°, 45cm length

Articulating atraumatic Diamond-Jaw® grasper, 5mm


	Catalog no.	Dimensions
	89-0509	40°, 34cm length
	89-0503	90°, 34cm length
	89-4509	40°, 45cm length
	89-4503	90°, 45cm length

Accessories


88-6281	Sterilization tray
88-6282	8-slot, 5mm fitting
88-6283	Single-slot, 5-10mm fitting
88-9861	Luer lock plug, 10/pack
89-6119	Protective sterilization sleeve, short
89-6124	Multi-positioning forceps tube




Articulating Diamond-Jaw grasper, 5mm

	Catalog no.	Dimensions
	89-0508	40°, 34cm length
	89-0502	90°, 34cm length
	89-4508	40°, 45cm length
	89-4502	90°, 45cm length

Articulating fenestrated grasper, 5mm

	Catalog no.	Dimensions
	89-0501	40°, 34cm length
	89-0505	90°, 34cm length
	89-4501	40°, 45cm length
	89-4505	90°, 45cm length

Articulating MARYLAND dissector, 5mm

	Catalog no.	Dimensions
	89-0511	40°, 34cm length
	89-0510	90°, 34cm length
	89-4511	40°, 45cm length
	89-4510	90°, 45cm length

Snowden-Pencer warranty

Snowden-Pencer articulating instruments feature a one-year service warranty, which covers the most common repairs laparoscopic instruments require, and a lifetime warranty against defects in materials and workmanship. Any product that proves to be defective in materials and workmanship will be repaired or replaced at our discretion. This guarantee (warranty) does not cover damage caused by abuse or misuse and is voided if personnel other than authorized Snowden-Pencer personnel or service resources repair or modify the instrument.

To place an order or for more information, please contact your
V. Mueller® and Snowden-Pencer Sales Representative, call Customer
Service at **800.323.9088** or visit carefusion.com/snowden-pencer

CareFusion
Waukegan, IL

carefusion.com

