

Wound drainage /
Surgery

SurgiLife

SURGILIFE Ε.Π.Ε.

ΝΟΣΟΚΟΜΕΙΑΚΑ ΧΕΙΡΟΥΡΓΙΚΑ ΥΛΙΚΑ

ΓΡΑΦΕΙΑ-ΛΟΓΙΣΤΗΡΙΟ: Σαρανταπόρου 15, ΑΘΗΝΑ Τ.Κ. 111 44


ΕΔΡΑ: Δροσίνι 2, ΑΘΗΝΑ Τ.Κ. 111 41

Τηλ. 210-21.16.291 - Fax. 210-21.16.292 - e-mail: surgilife@surgilife.gr

PRODUCT CATALOG



Primed[®]



SurgiLife

SURGILIFE E.Π.Ε.

ΝΟΣΟΚΟΜΕΙΑΚΑ ΧΕΙΡΟΥΡΓΙΚΑ ΥΛΙΚΑ

ΓΡΑΦΕΙΑ-ΛΟΓΙΣΤΗΡΙΟ: Σοφραντισπόρου 15, ΑΘΗΝΑ Τ.Κ. 111 44

ΕΔΡΑ: Δροσίνης 2, ΑΘΗΝΑ Τ.Κ. 111 41

Τηλ. 210-21.16.291 - Fax. 210-21.16.292 - e-mail: surgilife@surgilife.gr




PRIVAC® 600 ml complete set

- Redon bottle, connecting tube, Redon drain with trocar (guide needle)

Ch.	Drain length	Connector	PU/SU	REF
8	500 mm	Large lock	35/35	22400
10	800 mm	Large lock	35/35	22401
12	800 mm	Large lock	35/35	22402
14	800 mm	Large lock	35/35	22403
16	800 mm	Large lock	35/35	22404
18	800 mm	Large lock	35/35	22405


Ch.	Drain length	Connector	PU/SU	REF
8	500 mm	Luer lock	35/35	24750
10	800 mm	Luer lock	35/35	24751
12	800 mm	Luer lock	35/35	24752
14	800 mm	Luer lock	35/35	24753
16	800 mm	Luer lock	35/35	24754
18	800 mm	Luer lock	35/35	24755



PRIVAC® 600 ml OR system

- Redon bottle, connecting tube

Connector	PU/SU	REF
Large lock	35/35	24504
Luer lock	35/35	22407



PRIVAC® 600 ml ward system (replacement bottle)

- Redon bottle

Connector	PU/SU	REF
Large lock	35/35	24505
Luer lock	35/35	22406



Large lock



Luer lock

CONNECTION
OPTIONS

LOW-VACUUM SYSTEMS
PRI-LOW-VAC® SAVE
500 ml

SurgiLife


SURGILIFE Ε.Π.Ε.

ΝΟΣΟΚΟΜΕΙΑΚΑ ΧΕΙΡΟΥΡΓΙΚΑ ΥΛΙΚΑ

ΓΡΑΦΕΙΑ-ΛΟΓΙΣΤΗΡΙΟ: Σαρανταπόρου 15, ΑΘΗΝΑ Τ.Κ. 111 44

ΕΔΡΑ: Δροσίνη 2, ΑΘΗΝΑ Τ.Κ. 111 41


Τηλ. 210-21.16.291 - Fax. 210-21.16.292 - e-mail: surgilife@surgilife.gr



PRI-LOW-VAC® SAVE 500 ml complete set

- Redon bellows with connecting tube, Redon drain (PVC) with trocar


Ch.	Drain length	PU/SU	REF
10	800 mm	35/ 35	21981
12	800 mm	35/ 35	21987
14	800 mm	35/ 35	21988
16	800 mm	35/ 35	21989
18	800 mm	35/ 35	21980



PRI-LOW-VAC® SAVE 500 ml OR system

- Redon bellows with connecting tube

PU/SU	REF
35/ 35	21972



PRI-LOW-VAC® SAVE 500 ml ward system

(replacement bellows)


- Redon bellows

PU/SU	REF
35/ 35	21971

SURGILIFE Ε.Π.Ε.
ΝΟΣΟΚΟΜΕΙΑΚΑ ΧΕΙΡΟΥΡΓΙΚΑ ΥΛΙΚΑ
ΓΡΑΦΕΙΑ-ΛΟΓΙΣΤΗΡΙΟ: Σαρανταπόρου 15, ΑΘΗΝΑ Τ.Κ. 111 44
ΕΔΡΑ: Δροσίνι 2, ΑΘΗΝΑ Τ.Κ. 111 41
Τηλ. 210-21.16.291 - Fax. 210-21.16.292 - e-mail: surgilife@surgilife.gr

Passive Drainage (Gravity drainage)


Passive drainage is used for open or half-open wound drainage by gravity. This drainage system consists of a drainage bag with non-return valve, a connecting tube and a silicone soft drain.



Product features

- Drainage bag 500 ml
 - With non-return valve to avoid secretion backflow
 - With large-lumen tap
 - With scale
 - Patient information can be entered in pre-printed fields
- Silicone soft drain
 - Soft, medical grade silicone ensures good tissue compatibility
 - Rounded, atraumatic tip
 - 6 side openings
 - With radiopaque stripe

Ch.	Drain length	PU/SU	REF
8	1000 mm	25/50	22050
10	1000 mm	25/50	22051
12	1000 mm	25/50	21910
14	1000 mm	25/50	21898
15	1000 mm	25/50	21911
16	1000 mm	25/50	21845
18	1000 mm	25/50	21912
20	1000 mm	25/50	21888
21	1000 mm	25/50	21913
24	1000 mm	25/50	21914
26	1000 mm	25/50	21134
27	1000 mm	25/50	21915
30	1000 mm	25/50	21916
33	1000 mm	25/50	21917
36	1000 mm	25/50	21918
39	1000 mm	25/50	21135



Soft drain